



## SQLE rabbit pAb

货号	YP-Ab-18125
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	SQLE
基因名称	SQLE ERG1
蛋白名称	Squalene monooxygenase (EC 1.14.13.132) (Squalene epoxidase) (SE)
免疫原	Synthesized peptide derived from human SQLE
特异性	This antibody detects endogenous levels of SQLE at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	63kD
功能	Catalyzes the stereospecific oxidation of squalene to (S)-2,3-epoxysqualene, and is considered to be a rate-limiting enzyme in steroid biosynthesis.
细胞定位	Microsome membrane ; Peripheral membrane protein . Endoplasmic reticulum membrane ; Peripheral membrane protein .
组织表达	Detected in liver (at protein level).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images