



DUSP5 rabbit pAb

货号	YP-Ab-18126
同位型	IgG
应用	WB
种属	Human,Rat
靶点	DUSP5
基因名称	DUSP5 VH3
蛋白名称	Dual specificity protein phosphatase 5 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH3)
免疫原	Synthesized peptide derived from human DUSP5
特异性	This antibody detects endogenous levels of DUSP5 at Human,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	42kD
功能	Dual specificity protein phosphatase; active with phosphotyrosine, phosphoserine and phosphothreonine residues. The highest relative activity is toward ERK1.
细胞定位	Nucleus .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images