

UPB1 rabbit pAb

货号	YP-Ab-18128
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	UPB1
基因名称	UPB1 BUP1
蛋白名称	Beta-ureidopropionase (EC 3.5.1.6) (BUP-1) (Beta-alanine synthase) (N-carbamoyl-beta-alanine amidohydrolase)
免疫原	Synthesized peptide derived from human UPB1
特异性	This antibody detects endogenous levels of UPB1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	42kD
功能	Catalyzes a late step in pyrimidine degradation . Converts N-carbamoyl-beta-alanine (3-ureidopropanoate) into beta-alanine, ammonia and carbon dioxide . Likewise, converts N-carbamoyl-beta-aminoisobutyrate (3-ureidoisobutyrate) into beta-aminoisobutyrate, ammonia and carbon dioxide (Probable).
细胞定位	Cytoplasm.
组织表达	Detected in liver (at protein level).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

