

## TSPAN6 rabbit pAb

货号	YP-Ab-18132
同位型	IgG
应用	WB
种属	Human
靶点	TSPAN6
基因名称	TSPAN6 TM4SF6 UNQ767/PRO1560
蛋白名称	Tetraspanin-6 (Tspan-6) (A15 homolog) (Putative NF-kappa-B-activating protein 321) (T245 protein) (Tetraspanin TM4-D) (Transmembrane 4 superfamily member 6)
免疫原	Synthesized peptide derived from human TSPAN6
特异性	This antibody detects endogenous levels of TSPAN6 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	27kD
细胞定位	Membrane; Multi-pass membrane protein.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images