

# LRRC1 rabbit pAb

货号	YP-Ab-18139
同位型	IgG
应用	WB
种属	
靶点	LRRC1
基因名称	Human;Mouse
蛋白名称	Leucine-rich repeat-containing protein 1 (LANO adapter protein) (LAP and no PDZ protein)
免疫原	Synthesized peptide derived from human LRRC1
特异性	This antibody detects endogenous levels of LRRC1 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	58kD
细胞定位	Cytoplasm . Membrane ; Peripheral membrane protein . Localized at the basolateral side of epithelial cells.
组织表达	Expressed strongly in testis and placenta, followed by heart, lung, kidney, thyroid, trachea, colon, prostate and pancreas.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images