



PLAC8 rabbit pAb

| | |
|--------|---|
| 货号 | YP-Ab-18150 |
| 同位型 | IgG |
| 应用 | WB |
| 种属 | Human;Mouse |
| 靶点 | PLAC8 |
| 基因名称 | PLAC8 BM-004 |
| 蛋白名称 | Placenta-specific gene 8 protein (Protein C15) |
| 免疫原 | Synthesized peptide derived from human PLAC8 |
| 特异性 | This antibody detects endogenous levels of PLAC8 at Human, Mouse |
| 组成 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| 来源 | Polyclonal, Rabbit,IgG |
| 稀释 | WB 1:500-2000 |
| 纯化工艺 | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| 分子量 | 13kD |
| 组织表达 | Expressed at high levels in plasmacytoid dendritic cells. High expression in spleen, lymph nodes, peripheral blood leukocytes, and bone marrow, with lower expression in thymus, appendix, and fetal liver. |
| 浓度 | 1 mg/ml |
| 储存 | -15°C to -25°C/1 year(Do not lower than -25°C) |
| 有关注意事项 | Avoid repeated freezing and thawing! |
| 使用建议 | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |

Products Images