



ECH1 rabbit pAb

货号	YP-Ab-18157
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	ECH1
基因名称	ECH1
蛋白名称	Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial (EC 5.3.3.-)
免疫原	Synthesized peptide derived from human ECH1
特异性	This antibody detects endogenous levels of ECH1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	36kD
功能	Isomerization of 3-trans,5-cis-dienoyl-CoA to 2-trans,4-trans-dienoyl-CoA.
细胞定位	Mitochondrion . Peroxisome .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images