



CHCHD5 rabbit pAb

货号	YP-Ab-18158
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	CHCHD5
基因名称	CHCHD5 C2orf9
蛋白名称	Coiled-coil-helix-coiled-coil-helix domain-containing protein 5
免疫原	Synthesized peptide derived from human CHCHD5
特异性	This antibody detects endogenous levels of CHCHD5 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	12kD
细胞定位	Mitochondrion intermembrane space .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images