



## HECTD3 rabbit pAb

货号	YP-Ab-18166
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	HECTD3
基因名称	HECTD3
蛋白名称	E3 ubiquitin-protein ligase HECTD3 (EC 6.3.2.-) (HECT domain-containing protein 3)
免疫原	Synthesized peptide derived from human HECTD3
特异性	This antibody detects endogenous levels of HECTD3 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	95kD
功能	E3 ubiquitin ligases accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Mediates ubiquitination of TRIOBP and its subsequent proteasomal degradation, thus facilitating cell cycle progression by regulating the turn-over of TRIOBP. Mediates also ubiquitination of STX8 (By similarity).
细胞定位	Cytoplasm, perinuclear region .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

