



# CPSF4 rabbit pAb

货号	YP-Ab-18167
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	CPSF4
基因名称	CPSF4 CPSF30 NAR NEB1
蛋白名称	Cleavage and polyadenylation specificity factor subunit 4 (Cleavage and polyadenylation specificity factor 30 kDa subunit) (CPSF 30 kDa subunit) (NS1 effector domain-binding protein 1) (Neb-1) (No arc
免疫原	Synthesized peptide derived from human CPSF4
特异性	This antibody detects endogenous levels of CPSF4 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	30kD
功能	Component of the cleavage and polyadenylation specificity factor (CPSF) complex that play a key role in pre-mRNA 3'-end formation, recognizing the AAUAAA signal sequence and interacting with poly(A) polymerase and other factors to bring about cleavage and poly(A) addition. CPSF4 binds RNA polymers with a preference for poly(U).
细胞定位	Nucleus.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

