



## PCSK2 rabbit pAb

货号	YP-Ab-18168
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	PCSK2
基因名称	PCSK2 NEC2
蛋白名称	Neuroendocrine convertase 2 (NEC 2) (EC 3.4.21.94) (KEX2-like endoprotease 2) (Prohormone convertase 2) (Proprotein convertase 2) (PC2)
免疫原	Synthesized peptide derived from human PCSK2
特异性	This antibody detects endogenous levels of PCSK2 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	70kD
功能	Involved in the processing of hormone and other protein precursors at sites comprised of pairs of basic amino acid residues. Responsible for the release of glucagon from proglucagon in pancreatic A cells.
细胞定位	Cytoplasmic vesicle, secretory vesicle. Secreted . Localized in the secretion granules.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images