



OIP5 rabbit pAb

货号	YP-Ab-18170
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	OIP5
基因名称	OIP5 MIS18B
蛋白名称	Protein Mis18-beta (Cancer/testis antigen 86) (CT86) (Opa-interacting protein 5) (OIP-5)
免疫原	Synthesized peptide derived from human OIP5
特异性	This antibody detects endogenous levels of OIP5 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	25kD
功能	Required for recruitment of CENPA to centromeres and normal chromosome segregation during mitosis.
细胞定位	Nucleus . Chromosome . Chromosome, centromere . Associated with centromeres in interphase cells, from late anaphase to the G1 phase. Not detected on centromeres during earlier phases of mitosis. Associated with chromatin.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



优品生物