

COMMD8 rabbit pAb

货号	YP-Ab-18171
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	COMMD8
基因名称	COMMD8 MDS022
蛋白名称	COMM domain-containing protein 8
免疫原	Synthesized peptide derived from human COMMD8
特异性	This antibody detects endogenous levels of COMMD8 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	20kD
功能	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes . May down-regulate activation of NF-kappa-B .
细胞定位	Cytoplasm . Nucleus .
组织表达	Widely expressed with highest expression in thyroid.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images