



NOTUM rabbit pAb

货号	YP-Ab-18172
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	NOTUM
基因名称	NOTUM OK/SW-CL.30
蛋白名称	Protein notum homolog (EC 3.-.-.)
免疫原	Synthesized peptide derived from human NOTUM
特异性	This antibody detects endogenous levels of NOTUM at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	55kD
功能	Carboxylesterase that acts as a key negative regulator of the Wnt signaling pathway by specifically mediating depalmitoleoylation of WNT proteins. Serine palmitoleoylation of WNT proteins is required for efficient binding to frizzled receptors .
细胞定位	Secreted .
组织表达	Rarely expressed in adult normal tissues.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

