



MSL3 rabbit pAb

货号	YP-Ab-18183
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	MSL3
基因名称	MSL3 MSL3L1
蛋白名称	Male-specific lethal 3 homolog (Male-specific lethal-3 homolog 1) (Male-specific lethal-3 protein-like 1) (MSL3-like 1)
免疫原	Synthesized peptide derived from human MSL3
特异性	This antibody detects endogenous levels of MSL3 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	57kD
功能	Has a role in chromatin remodeling and transcriptional regulation . Has a role in X inactivation . Component of the MSL complex which is responsible for the majority of histone H4 acetylation at 'Lys-16' which is implicated in the formation of higher-order chromatin structure . Specifically recognizes histone H4 monomethylated at 'Lys-20' (H4K20Me1) in a DNA-dependent manner and is proposed to be involved in chromosomal targeting of the MSL complex .
细胞定位	Nucleus .
组织表达	Expressed in many tissues including liver, pancreas, heart, lung, kidney, skeletal muscle, brain, and placenta, with highest expression in skeletal muscle and heart.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

