

# CTHRC1 rabbit pAb

货号	YP-Ab-18199
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	CTHRC1
基因名称	CTHRC1 UNQ762/PRO1550
蛋白名称	Collagen triple helix repeat-containing protein 1 (Protein NMTC1)
免疫原	Synthesized peptide derived from human CTHRC1
特异性	This antibody detects endogenous levels of CTHRC1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	27kD
功能	May act as a negative regulator of collagen matrix deposition.
细胞定位	Secreted, extracellular space, extracellular matrix .
组织表达	Isoform 1 is expressed in calcified atherosclerotic plaque and chondrocyte-like cells.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images