



STXBP3 rabbit pAb

货号	YP-Ab-18218
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	STXBP3
基因名称	STXBP3
蛋白名称	Syntaxin-binding protein 3 (Platelet Sec1 protein) (PSP) (Protein unc-18 homolog 3) (Unc18-3) (Protein unc-18 homolog C) (Unc-18C)
免疫原	Synthesized peptide derived from human STXBP3
特异性	This antibody detects endogenous levels of STXBP3 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	65kD
功能	Together with STX4 and VAMP2, may play a role in insulin-dependent movement of GLUT4 and in docking/fusion of intracellular GLUT4-containing vesicles with the cell surface in adipocytes.
细胞定位	Cytoplasm, cytosol . Cell membrane . In platelets, predominantly cytosolic. Low amounts membrane-associated.
组织表达	Megakaryocytes and platelets.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

