



TET1 Rabbit pAb

Catalog No	YP-Ab-18231
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC-F, IHC-P, ICC/IF, FC, ELISA
Gene Name	TET
Alternative Names	Leukemia associated protein with a CXXC domain; CXXC 6; CXXC finger 6; CXXC type zinc finger protein 6; CXXC-type zinc finger protein 6; CXXC6; KIAA1676; LCX; Leukemia-associated protein with a CXXC domain; Methylcytosine dioxygenase TET1; Ten eleven translocation 1 gene protein; Ten eleven translocation 1 gene protein homolog; Ten-eleven translocation 1 gene protein; Tet 1; Tet oncogene 1; TET1; TET1_HUMAN; TET1/CXXC6.
Research Field	
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 235kDa
Clonality	Polyclonal
Clonality No.	
Dilution	IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/5000-1/10000
Immunogen	KLH conjugated synthetic peptide derived from human TET1
Purification	affinity purified by Protein A
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images