



SUMF2 rabbit pAb

货号	YP-Ab-18234
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	SUMF2
基因名称	SUMF2 PSEC0171 UNQ1968/PRO4500
蛋白名称	Sulfatase-modifying factor 2 (C-alpha-formylglycine-generating enzyme 2)
免疫原	Synthesized peptide derived from human SUMF2
特异性	This antibody detects endogenous levels of SUMF2 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	33kD
功能	Lacks formylglycine generating activity and is unable to convert newly synthesized inactive sulfatases to their active form. Inhibits the activation of sulfatases by SUMF1.
细胞定位	Endoplasmic reticulum lumen .
组织表达	Detected in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Highest levels in kidney, liver and placenta.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

