



## FLJ14213 rabbit pAb

货号	YP-Ab-18241
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	PRR5L
基因名称	PRR5L PROTOR2
蛋白名称	Proline-rich protein 5-like (Protein observed with Rictor-2) (Protor-2)
免疫原	Synthesized peptide derived from human FLJ14213
特异性	This antibody detects endogenous levels of FLJ14213 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	40kD
功能	Associates with the mTORC2 complex that regulates cellular processes including survival and organization of the cytoskeleton . Regulates the activity of the mTORC2 complex in a substrate-specific manner preventing for instance the specific phosphorylation of PKCs and thereby controlling cell migration . Plays a role in the stimulation of ZFP36-mediated mRNA decay of several ZFP36-associated mRNAs, such as TNF-alpha and GM-CSF, in response to stress . Required for ZFP36 localization to cytoplasmic stress granule (SG) and P-body (PB) in response to stress .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

