



## SASS6 rabbit pAb

货号	YP-Ab-18242
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	SASS6
基因名称	SASS6 SAS6
蛋白名称	Spindle assembly abnormal protein 6 homolog (HsSAS-6)
免疫原	Synthesized peptide derived from human SASS6
特异性	This antibody detects endogenous levels of SASS6 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	72kD
功能	Central scaffolding component of the centrioles ensuring their 9-fold symmetry. Required for centrosome biogenesis and duplication: required both for mother-centriole-dependent centriole duplication and deuterosome-dependent centriole amplification in multiciliated cells. Overexpression results in excess foci-bearing centriolar markers. Required for the recruitment of STIL to the procentriole and for STIL-mediated centriole amplification .
细胞定位	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Component of the deuterosome, a structure that promotes de novo centriole amplification in multiciliated cells that can generate more than 100 centrioles (By similarity). Component of the centrosome. Associated only transiently with nascent procentrioles during centriole biogenesis. .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

