



FAAH rabbit pAb

货号	YP-Ab-18243
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	FAAH
基因名称	FAAH FAAH1
蛋白名称	Fatty-acid amide hydrolase 1 (EC 3.5.1.99) (Anandamide amidohydrolase 1) (Oleamide hydrolase 1)
免疫原	Synthesized peptide derived from human FAAH
特异性	This antibody detects endogenous levels of FAAH at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	64kD
功能	Catalyzes the hydrolysis of endogenous amidated lipids like the sleep-inducing lipid oleamide ((9Z)-octadecenamide), the endocannabinoid anandamide (N-(5Z,8Z,11Z,14Z-eicosatetraenoyl)-ethanolamine), as well as other fatty amides, to their corresponding fatty acids, thereby regulating the signaling functions of these molecules . Hydrolyzes polyunsaturated substrate anandamide preferentially as compared to monounsaturated substrates . It can also catalyze the hydrolysis of the endocannabinoid 2-arachidonoylglycerol (2-(5Z,8Z,11Z,14Z-eicosatetraenoyl)-glycerol) . FAAH cooperates with PM20D1 in the hydrolysis of amino acid-conjugated fatty acids such as N-fatty acyl glycine and N-fatty acyl-L-serine, thereby acting as a physiological regulator of specific subsets of intracellular, but not of extracellular, N-fatty acyl amino acids (By similarity).
细胞定位	Endomembrane system ; Single-pass membrane protein . Cytoplasm, cytoskeleton . Seems to be attached to intracellular membranes and a portion of the cytoskeletal network.
组织表达	Highly expressed in the brain, small intestine, pancreas, skeletal muscle and testis. Also expressed in the kidney, liver, lung, placenta and prostate.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!

使用建议

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images