



ATL3 rabbit pAb

货号	YP-Ab-18254
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	ATL3
基因名称	ATL3
蛋白名称	Atlastin-3 (EC 3.6.5.-)
免疫原	Synthesized peptide derived from human ATL3
特异性	This antibody detects endogenous levels of ATL3 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	60kD
功能	GTPase tethering membranes through formation of trans-homooligomers and mediating homotypic fusion of endoplasmic reticulum membranes. Functions in endoplasmic reticulum tubular network biogenesis .
细胞定位	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Localizes to endoplasmic reticulum tubules and accumulates in punctuate structures corresponding to 3-way junctions, which represent crossing-points at which the tubules build a polygonal network. .
组织表达	Expressed in the central nervous system and in dorsal root ganglia neurons. Expressed in peripheral tissues (at protein level).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

