



SNRNP200 rabbit pAb

货号	YP-Ab-18257
同位型	IgG
应用	WB
种属	Human
靶点	SNRNP200
基因名称	SNRNP200 ASCC3L1 HELIC2 KIAA0788
蛋白名称	U5 small nuclear ribonucleoprotein 200 kDa helicase (EC 3.6.4.13) (Activating signal cointegrator 1 complex subunit 3-like 1) (BRR2 homolog) (U5 snRNP-specific 200 kDa protein) (U5-200KD)
免疫原	Synthesized peptide derived from human SNRNP200
特异性	This antibody detects endogenous levels of SNRNP200 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	235kD
功能	Plays role in pre-mRNA splicing as core component of precatalytic, catalytic and postcatalytic spliceosomal complexes . Involved in spliceosome assembly, activation and disassembly. Mediates changes in the dynamic network of RNA-RNA interactions in the spliceosome. Catalyzes the ATP-dependent unwinding of U4/U6 RNA duplexes, an essential step in the assembly of a catalytically active spliceosome.
细胞定位	Nucleus .
组织表达	Widely expressed.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

