



DIMT1L rabbit pAb

货号	YP-Ab-18263
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	DIMT1
基因名称	DIMT1 DIMT1L HUSSY-05
蛋白名称	Probable dimethyladenosine transferase (EC 2.1.1.183) (DIM1 dimethyladenosine transferase 1 homolog) (DIM1 dimethyladenosine transferase 1-like) (Probable 18S rRNA (adenine(1779)-N(6)/adenine(1780)-N(
免疫原	Synthesized peptide derived from human DIMT1L
特异性	This antibody detects endogenous levels of DIMT1L at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	34kD
功能	Specifically dimethylates two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18S rRNA in the 40S particle . Involved in the pre-rRNA processing steps leading to small-subunit rRNA production independently of its RNA-modifying catalytic activity .
细胞定位	Nucleus, nucleoplasm . Nucleus, nucleolus .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

