

PNO1 rabbit pAb

货号	YP-Ab-18264
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	PNO1
基因名称	PNO1
蛋白名称	RNA-binding protein PNO1
免疫原	Synthesized peptide derived from human PNO1
特异性	This antibody detects endogenous levels of PNO1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	28kD
功能	Positively regulates dimethylation of two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18S rRNA .
细胞定位	Nucleus, nucleolus .
组织表达	Expressed in liver, lung, spleen and kidney. Weakly expressed in thymus, testis and ovary. Weakly or not expressed in heart, brain, skeletal muscle, placenta, pancreas, prostate, small intestine, colon and peripheral blood leukocytes.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

