



NOB1 rabbit pAb

货号	YP-Ab-18265
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	NOB1
基因名称	NOB1 ART4 NOB1P PSMD8BP1 MSTP158
蛋白名称	RNA-binding protein NOB1 (Phosphorylation regulatory protein HP-10) (Protein ART-4)
免疫原	Synthesized peptide derived from human NOB1
特异性	This antibody detects endogenous levels of NOB1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	45kD
功能	May play a role in mRNA degradation (Probable). Endonuclease required for processing of 20S pre-rRNA precursor and biogenesis of 40S ribosomal subunits (By similarity).
细胞定位	Nucleus .
组织表达	Detected in liver, lung, placenta, endothelial cells and spleen.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

