



## CILP2 rabbit pAb

货号	YP-Ab-18273
同位型	IgG
应用	WB
种属	Human
靶点	CILP2
基因名称	CILP2
蛋白名称	Cartilage intermediate layer protein 2 (CILP-2) [Cleaved into: Cartilage intermediate layer protein 2 C1; Cartilage intermediate layer protein 2 C2]
免疫原	Synthesized peptide derived from human CILP2
特异性	This antibody detects endogenous levels of CILP2 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	127kD
功能	May play a role in cartilage scaffolding.
细胞定位	Secreted, extracellular space, extracellular matrix .
组织表达	Expressed in articular chondrocytes but not in knee meniscal cartilage cells (PubMed:12746903). Localizes to the intermediate to deep zone of articular cartilage (PubMed:21880736).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

