

MSRA rabbit pAb

| | |
|--------|--|
| 货号 | YP-Ab-18285 |
| 同位型 | IgG |
| 应用 | WB |
| 种属 | Human;Mouse;Rat |
| 靶点 | MSRA |
| 基因名称 | MSRA |
| 蛋白名称 | Mitochondrial peptide methionine sulfoxide reductase (EC 1.8.4.11) (Peptide-methionine (S)-S-oxide reductase) (Peptide Met(O) reductase) (Protein-methionine-S-oxide reductase) (PMSR) |
| 免疫原 | Synthesized peptide derived from human MSRA |
| 特异性 | This antibody detects endogenous levels of MSRA at Human, Mouse,Rat |
| 组成 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| 来源 | Polyclonal, Rabbit,IgG |
| 稀释 | WB 1:500-2000 |
| 纯化工艺 | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| 分子量 | 26kD |
| 功能 | Has an important function as a repair enzyme for proteins that have been inactivated by oxidation. Catalyzes the reversible oxidation-reduction of methionine sulfoxide in proteins to methionine. |
| 细胞定位 | [Isoform 1]: Mitochondrion.; [Isoform 2]: Cytoplasm.; [Isoform 3]: Cytoplasm. Nucleus.; [Isoform 5]: Cytoplasm. Membrane ; Lipid-anchor . |
| 组织表达 | Ubiquitous. Highest expression in adult kidney and cerebellum, followed by liver, heart ventricles, bone marrow and hippocampus. |
| 浓度 | 1 mg/ml |
| 储存 | -15°C to -25°C/1 year(Do not lower than -25°C) |
| 有关注意事项 | Avoid repeated freezing and thawing! |
| 使用建议 | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images

