



## KIFC1 rabbit pAb

货号	YP-Ab-18292
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	KIFC1
基因名称	KIFC1 HSET KNSL2
蛋白名称	Kinesin-like protein KIFC1 (Kinesin-like protein 2) (Kinesin-related protein HSET)
免疫原	Synthesized peptide derived from human KIFC1
特异性	This antibody detects endogenous levels of KIFC1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	74kD
功能	Minus end-directed microtubule-dependent motor required for bipolar spindle formation . May contribute to movement of early endocytic vesicles (By similarity). Regulates cilium formation and structure (By similarity).
细胞定位	Nucleus . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle . Early endosome . Associated with nucleus during interphase, centrosomes in early and spindle in later mitosis. .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images