



KIAA1522蛋白抗体

货号	YP-Ab-18309
英文名称	Anti-KIAA1522
别名	K1522 HUMAN; Kiaa1522;Uncharacterized proteinKAA1522.
交叉反应	Human,Mouse, Rat,Chicken, Dog,Pig,Cow, Horse,Rabbit
性状	Lyophilized or Liquid
免疫原	KLHconjugated synthetic peptide derived fromhuman KIAA1522
来源	Rabbit
克隆类型	polyclonal
分子量	114kDa
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images