




ERM Monoclonal Antibody

Catalog No	YP-Ab-00971
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	ETV5
Protein Name	ETS translocation variant 5
Immunogen	Purified recombinant fragment of human ERM expressed in E. Coli.
Specificity	ERM Monoclonal Antibody detects endogenous levels of ERM protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ETV5; ERM; ETS translocation variant 5; Ets-related protein ERM
Observed Band	
Cell Pathway	Nucleus.
Tissue Specificity	Ubiquitous.
Function	function: Binds to DNA sequences containing the consensus nucleotide core sequence GGAA., mass spectrometry: PubMed:11840567, similarity: Belongs to the ETS family., similarity: Contains 1 ETS DNA-binding domain., tissue specificity: Ubiquitous.,
Background	function: Binds to DNA sequences containing the consensus nucleotide core sequence GGAA., mass spectrometry: PubMed:11840567, similarity: Belongs to the ETS family., similarity: Contains 1 ETS DNA-binding domain., tissue specificity: Ubiquitous.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using ERM Monoclonal Antibody against Jurkat (1), NIH/3T3 (2) and MCF-7 (3) cell lysate.