




DQX1 Polyclonal Antibody


Catalog No	YP-Ab-01669
Isotype	IgG
Reactivity	Human;Mouse;Monkey
Applications	WB;ELISA
Gene Name	DQX1
Protein Name	ATP-dependent RNA helicase DQX1
Immunogen	The antiserum was produced against synthesized peptide derived from human DQX1. AA range:551-600
Specificity	DQX1 Polyclonal Antibody detects endogenous levels of DQX1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1
Observed Band	79kD
Cell Pathway	Nucleus .
Tissue Specificity	Placenta,
Function	similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,
Background	similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DQX1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).



Western blot analysis of lysates from COS7 cells, using DQX1 Antibody. The lane on the right is blocked with the synthesized peptide.