





Prolactin mouse mAb(ABT220)

Catalog No	YP-Ab-15617
Isotype	IgG
Reactivity	Human
Applications	IHC;IF
Gene Name	PRL
Protein Name	Prolactin
Immunogen	Synthesized peptide derived from human Prolactin
Specificity	The antibody can specifically recognize human Prolactin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.19% sodium azide.
Source	Mouse, Monoclonal/IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Prolactin (PRL)
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	Cytoplasmic
Function	function:Prolactin acts primarily on the mammary gland by promoting lactation.,online information:Prolactin entry,similarity:Belongs to the somatotropin/prolactin family.,
Background	This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Aug 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


Products Images



Human pituitary adenoma tissue was stained with Anti-Prolactin (ABT220) Antibody



Human pituitary adenoma tissue was stained with Anti-Prolactin (ABT220) Antibody



Human pituitary adenoma tissue was stained with Anti-Prolactin (ABT220) Antibody