




Thyroperoxidase rabbit pAb

Catalog No	YP-Ab-15720
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC
Gene Name	TPO
Protein Name	Thyroperoxidase
Immunogen	Synthesized peptide derived from human Thyroperoxidase AA range: 651-700
Specificity	This antibody detects endogenous levels of Human Thyroperoxidase
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Thyroid peroxidase (TPO;EC 1.11.1.8)
Observed Band	
Cell Pathway	Membrane; Single-pass type I membrane protein.; [Isoform 3]: Cell surface.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).