



EFNA3 (Cleaved-Gly214) rabbit pAb

Catalog No	YP-Ab-15820
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	EFNA3 EFL2 EPLG3 LERK3
Protein Name	EFNA3 (Cleaved-Gly214)
Immunogen	Synthesized peptide derived from human EFNA3 (Cleaved-Gly214)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat EFNA3 (Cleaved-Gly214, protein was cleaved amino acid sequence between 214-215)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ephrin-A3 (EFL-2;EHK1 ligand;EHK1-L;EPH-related receptor tyrosine kinase ligand 3;LERK-3)
Observed Band	24kD
Cell Pathway	Cell membrane; Lipid-anchor, GPI-anchor.
Tissue Specificity	Expressed in brain, skeletal muscle, spleen, thymus, prostate, testis, ovary, small intestine, and peripheral blood leukocytes.
Function	cell-cell signaling,
Background	similarity:Belongs to the ephrin family.,tissue specificity:Expressed in brain, skeletal muscle, spleen, thymus, prostate, testis, ovary, small intestine, and peripheral blood leukocytes.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images