




# CC157 rabbit pAb

|                                  |                                                                                                                                                         |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-Ab-12281                                                                                                                                             |
| <b>Isotype</b>                   | IgG                                                                                                                                                     |
| <b>Reactivity</b>                | Human; Mouse                                                                                                                                            |
| <b>Applications</b>              | WB                                                                                                                                                      |
| <b>Gene Name</b>                 | CCDC157 KIAA1656                                                                                                                                        |
| <b>Protein Name</b>              | CC157                                                                                                                                                   |
| <b>Immunogen</b>                 | Synthesized peptide derived from human CC157 AA range: 16-66                                                                                            |
| <b>Specificity</b>               | This antibody detects endogenous levels of CC157 at Human/Mouse                                                                                         |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                 |
| <b>Source</b>                    | Polyclonal, Rabbit,IgG                                                                                                                                  |
| <b>Purification</b>              | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                               |
| <b>Dilution</b>                  | WB 1: 500-2000                                                                                                                                          |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                                 |
| <b>Purity</b>                    | ≥90%                                                                                                                                                    |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                            |
| <b>Synonyms</b>                  |                                                                                                                                                         |
| <b>Observed Band</b>             |                                                                                                                                                         |
| <b>Cell Pathway</b>              |                                                                                                                                                         |
| <b>Tissue Specificity</b>        |                                                                                                                                                         |
| <b>Function</b>                  | sequence caution:Intron retention. This sequence is incomplete at the 5'- and 3'-ends and extensively differs from that shown due to multiple introns., |
| <b>Background</b>                |                                                                                                                                                         |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                                    |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                       |



## Products Images



Western blot analysis of lysates from U2OS cells,  
primary antibody was diluted at 1:1000, 4° over night