




RNK rabbit pAb

Catalog No	YP-Ab-12328
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	RNASEK
Protein Name	RNK
Immunogen	Synthesized peptide derived from human RNK AA range: 27-77
Specificity	This antibody detects endogenous levels of RNK at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1: 50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endomembrane system ; Multi-pass membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed.
Function	function:Endoribonuclease which preferentially cleaves ApU and ApG phosphodiester bonds. Hydrolyzes UpU bonds at a lower rate.,similarity:Belongs to the RNase K family.,tissue specificity:Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).