




CO029 rabbit pAb

Catalog No	YP-Ab-12349
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C15orf29
Protein Name	CO029
Immunogen	Synthesized peptide derived from human CO029 AA range: 243-293
Specificity	This antibody detects endogenous levels of CO029 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cytoplasm, cytoskeleton, spindle pole . Localizes to the spindle poles only during mitosis. Sequestered to the nucleus during interphase. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night