




MPND rabbit pAb

Catalog No	YP-Ab-12353
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MPND
Protein Name	MPND
Immunogen	Synthesized peptide derived from human MPND AA range: 353-403
Specificity	This antibody detects endogenous levels of MPND at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	domain:The JAMM motif may mediate the protease activity.,function:Probable protease.,similarity:Belongs to the peptidase M67 family.,similarity:Contains 1 MPN (JAB/Mov34) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night