




NCKP5 rabbit pAb

Catalog No	YP-Ab-12405
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NCKAP5 ERIH NAP5
Protein Name	NCKP5
Immunogen	Synthesized peptide derived from human NCKP5 AA range: 965-1015
Specificity	This antibody detects endogenous levels of NCKP5 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Expressed in fetal and adult brain, leukocytes and fetal fibroblasts.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night