



glrx5 rabbit pAb

Catalog No	YP-Ab-12584
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GLRX5 C14orf87
Protein Name	glrx5
Immunogen	Synthesized peptide derived from human glrx5
Specificity	This antibody detects endogenous levels of glrx5 at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.76% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glutaredoxin-related protein 5, mitochondrial (Monothiol glutaredoxin-5)
Observed Band	
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

