




FKBPL Polyclonal Antibody


Catalog No	YP-Ab-13928
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	FKBPL
Protein Name	FK506-binding protein-like
Immunogen	The antiserum was produced against synthesized peptide derived from human FKBPL. AA range:251-300
Specificity	FKBPL Polyclonal Antibody detects endogenous levels of FKBPL protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FKBPL; DIR1; NG7; FK506-binding protein-like; WAF-1/CIP1 stabilizing protein 39; WISp39
Observed Band	38kD
Cell Pathway	endoplasmic reticulum membrane,cytosol,
Tissue Specificity	Ubiquitously expressed with higher levels in testis.
Function	function:May be involved in response to X-ray. Regulates p21 protein stability by binding to Hsp90 and p21.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 3 TPR repeats.,subunit:Forms a ternary complex with CDKN1A/p21 and HSP90AB1/Hsp90.,tissue specificity:Ubiquitously expressed with higher levels in testis.,
Background	The protein encoded by this gene has similarity to the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. The encoded protein is thought to have a potential role in the induced radioresistance. Also it appears to have some involvement in the control of the cell cycle. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Jurkat cells, using FKBPL Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).