



c-Myc Monoclonal Antibody

Catalog No	YP-Ab-00961
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	MYC
Protein Name	Myc proto-oncogene protein
Immunogen	Purified recombinant fragment of human c-Myc expressed in E. Coli.
Specificity	c-Myc Monoclonal Antibody detects endogenous levels of c-Myc protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MYC; BHLHE39; Myc proto-oncogene protein; Class E basic helix-loop-helix protein 39; bHLHe39; Proto-oncogene c-Myc; Transcription factor p64
Observed Band	
Cell Pathway	Nucleus, nucleoplasm . Nucleus, nucleolus .
Tissue Specificity	Cervix, Epithelium, Leukemia, Placenta, Promyelocytic I
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

