



p47 phox (Phospho Ser345) rabbit pAb

Catalog No	YP-Ab-02438
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	NCF1 NOXO2 SH3PXD1A
Protein Name	p47 phox (Phospho Ser345)
Immunogen	Synthesized peptide derived from human p47 phox (Phospho Ser345)
Specificity	This antibody detects endogenous levels of Human p47 phox (Phospho Ser345)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Neutrophil cytosol factor 1 (NCF-1;47 kDa autosomal chronic granulomatous disease protein;47 kDa neutrophil oxidase factor;NCF-47K;Neutrophil NADPH oxidase factor 1;Nox organizer 2;Nox-organizing protein 2;SH3 and PX domain-containing protein 1A;p47-phox)
Observed Band	69
Cell Pathway	Cytoplasm, cytosol . Membrane ; Peripheral membrane protein ; Cytoplasmic side
Tissue Specificity	Detected in peripheral blood monocytes and neutrophils (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images