



# PDK1 (Phospho Tyr243) rabbit pAb

<b>Catalog No</b>	YP-Ab-02440
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB;ELISA;IHC-p
<b>Gene Name</b>	PDK1 PDHK1
<b>Protein Name</b>	PDHK1 (Phospho Tyr243)
<b>Immunogen</b>	Synthesized peptide derived from human PDHK1 (Phospho Tyr243)
<b>Specificity</b>	This antibody detects endogenous levels of Human,Mouse,Rat PDHK1 (Phospho Tyr243)
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	[Pyruvate dehydrogenase [lipoamide]] kinase isozyme 1, mitochondrial (EC 2.7.11.2;Pyruvate dehydrogenase kinase isoform 1)
<b>Observed Band</b>	50
<b>Cell Pathway</b>	Mitochondrion matrix .
<b>Tissue Specificity</b>	Expressed predominantly in the heart. Detected at lower levels in liver, skeletal muscle and pancreas.
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

