



His-tag mouse mAb

Catalog No	YP-Ab-04505
Isotype	IgG
Reactivity	
Applications	WB
Gene Name	His-Tag
Protein Name	His Tag
Immunogen	Synthetic peptide corresponding to fragment (HHHHHH) His Tag
Specificity	This antibody specifically recognizes the 6X His-Tag at C-terminal or N-terminal of recombinant proteins, which can be used to detect His-Tag fusion protein of various expression vector.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	wb dilution 1:1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	6 His epitope tag; Hexa His tag; HHHHHH epitope tag; HHHHHH tag; His tag; Polyhistidine Tag.
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images