



Flag-Tag Monoclonal Antibody(2C5), AbFluor™ 488 Conjugated

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-04546 |
| Isotype | IgG |
| Reactivity | Species independent |
| Applications | WB;IP;IF/ICC |
| Gene Name | Flag Tag |
| Protein Name | Flag-Tag DDDDK TAG, DDDDK-TAG |
| Immunogen | |
| Specificity | Flag-Tag Monoclonal Antibody(2C5) AbFluor™ 488 Conjugated specially designed for your Immunofluorescence analysis. |
| Formulation | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. |
| Source | Monoclonal, Mouse IgG |
| Purification | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Dilution | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:2000. |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | Flag |
| Observed Band | |
| Cell Pathway | |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images