



K2C1B Polyclonal Antibody

Catalog No	YP-Ab-06645
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	KRT77 KRT1B
Protein Name	Keratin, type II cytoskeletal 1b (Cytokeratin-1B) (CK-1B) (Keratin-77) (K77) (Type-II keratin Kb39)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	K2C1B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	63
Cell Pathway	cytoskeleton,intermediate filament,keratin filament,extracellular exosome,
Tissue Specificity	Expressed exclusively in skin.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

