



TIG1 Polyclonal Antibody

Catalog No	YP-Ab-07723
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	RARRES1 TIG1
Protein Name	Retinoic acid receptor responder protein 1 (RAR-responsive protein TIG1) (Tazarotene-induced gene 1 protein)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TIG1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	32
Cell Pathway	Membrane ; Single-pass type III membrane protein .
Tissue Specificity	Skin,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

